


ARTIFICIAL LARYNX



INHEALTH[®]
TECHNOLOGIES

Sale & Service by:

N.R. Enterprises

F-11A (G-8 AREA) NEAR BSES OFFICE
HARI NAGAR, NEW DELHI-110 064, INDIA
PH.: 47027049, 47000230, 40079990, 9810155141
pradeep6goel@yahoo.com • www.nrenterprisesindia.com

voice restoration systems

An electrolarynx is a small device used to help you talk after a total laryngectomy. This device is hand held and powered by batteries.